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Thank You for buying this UNIT product!

Attention! When you getting a UNIT product you must obtain from the seller completed guarantee card. Only if you have it you can use the manufacturer`s warranty in all UNIT service centers.

READ ALL THE INSTRUCTIONS BEFORE USING.

PURPOSE OF THE DEVICE

The device is intended for measurement of the weight and of the percentage of body fat and hydration.

INSTRUCTION AND PRECAUTIONS

1. Using this device you can determine the percentage of body fat as well as the hydration. The principle of operation of these scales is based on Bio-impedance analysis (Bioelectric impedance analysis – BIA), which allows you to determine the percentage of body fat as well as the content of water in your organism based at such parameters as age, gender, height.
2. The data obtained during the measurements can be used only as reference information. It is not recommended to use these values for medical goals. For getting more information on the recommended diet, etc. see you doctor.
3. For people using pacemakers or any other implants with electronic components it is not recommended to use this device. During the using of the device the implant can be disabled.
4. To avoid accidental interferences do not use the scales near strong electromagnetic field. Always install and use the scales away from such devices.
5. Use a soft and dry cloth to clean the scales. Do not immerse the appliance in water. Do not use chemical cleaners to clean the appliance.
6. Before making measurements make sure that the working surface of the scale is dry.
7. Place the scales on the dry flat surface (not on the carpet). If you do not plan to use the devise for long time remove the batteries from it.

8. Do not drop the appliance and do not jump on it. Protect the scales from strokes, dust, chemical substances and sharp changes in temperature.
9. Do not insert foreign objects into the openings of the device.
10. Repair of scales can be carried out only in the authorized UNIT service centers.
11. If the device does not turn on, replace the batteries. If after installing new batteries the appliance does not switch on, contact an authorized service center.
12. If the error is on the display for long time, remove the batteries from the scales for 3 seconds. Then put them back and start the measurement again.
13. In some cases, data from the memory cells of your scales may be lost. It is not a warranty case as the data can be deleted from the memory of the scales due to the influence of electromagnetic waves from various electrical appliances.
14. The manufacturer is not responsible for non-compliance of the precautions, which causes damage to the health of users of the scales or third parties.

FEATURES

PARAMETER	DEFAULT SETTING	RANGE
Memory code	1	1 to 12
Gender	Male	Male / Female
Height	165 cm	100 cm – 250 cm
Age	25	10 – 100

Max. Weighing Capacity: 150 kg

Min. Weighing at: 2.0 kg

Weight Measuring Division: 100 g

Body Fat Measure Range: 4.0-60.0%

Body Fat Differential Resolution in %: 0.1%

Hydration Measure Range: 27.5-66.0%

Hydration Differential Resolution in %: 0.1%

Muscle Range: 20.0-56.0%

Muscle Division: 0.1%

12 Persons Memory

Auto On


Low Battery Power Indication ('LO') / Overload Indication ('Err') Power: 2 x AAA battery

LCD-DISPLAY DESCRIPTION


INITIAL PREPARATION

- SET Turn on the scale and prompt to next parameter. (Please refer to the Setting Personal Parameter overleaf).
- ▲ Advance the value by one interval each press. Continuous pressing will speed up the advance of value. For choosing the code, gender, age and height.
 - ▼ Decrease the value by one interval each press. Continuous pressing will speed up the Decrease of value. For choosing the code, gender, age and height.
1. Install battery into the compartment underneath the main scales unit.
 2. Place the unit on a dry, flat and solid ground surface, ensure the unit is standing stable. Avoid carpet, or wavy ground surface.

PARAMETER SETTING

OPERATING THE APPLIANCE

MEASUREMENT OF BODY FAT, MUSCLE AND BONE MASS PERCENTAGE

1. Turn on the scales. The LCD-display will light and data from the memory will be displayed. By pressing the «SET» button you can step-by-step set the required parameters for obtain the result:
 - a) Press the button (▼▲) to select the memory location where you want to store you parameters or which you are already using. Press «SET» again for confirmation of choosing parameter. The next parameter of installation will be displayed (Gender).
 - b) Press the button (▼▲) to select your gender. Press «SET» again for confirmation of choosing parameter. The next parameter of installation will be displayed (Height).
 - c) Press the button (▼▲) to select your height. Press «SET» again for confirmation of choosing parameter. The next parameter of installation will be displayed (Age).
 - d) Press the button (▼▲) to select your age. Press «SET» again for confirmation of choosing parameter. The settings have been completed!

You can press the «SET» button a few more times to change parameters you put down incorrectly (the cycle described above will b repeated again).

2. Stand on the working surface of the scales. You feet should be bare and dry. Make sure your feel are located on a metal surface (each one on two squares), make sure your feet do not touch each other (in order to avoid short circuit).

If you step on the surface of the scales and the LSD-display automatically turned off, press the «SET» button and wait until current settings will appear on the display.

3. As soon as you step on the scales, the display will show the weight indicator (first time it will be increasing until it fixes on the exact value).
4. Within a few seconds, the scales will automatically determine the percentage of your body fat, hydration, muscle and bone mass and also daily cal-

orie recommendation. Then the scales will show these results 5 times and automatically shut down.

During the determination of your body fat percentage you will see an analysis of the measurement results in the form of series of vertical symbols. This is given for indication of your Fat percentage level, telling where your fat level standing is comparing to the statistics norms. Whether you are considered as underweight, normal, or overweight. You can take this indication as a reference data of your condition.

NOTE: *Each time you press the «SET» button you can start the measurement again (including changing the settings of the previous user) in accordance with the above points.*

WEIGHT MEASUREMENT

5. To activate the scales click on the segment of the working surface. After a few seconds the value «0,0 kg» will appear.
6. Stand on the working surface gently (do not jump on the scales). As soon as you step on the scales there will be the weight indicator on the display (first time it will be increasing until it fixes on the exact value). After a few seconds the scales automatically turns off.

MEASUREMENT ERRORS

If there are errors during operation the display shows «Err». In this case wait until the scales is turned off to start the measurement again, whether look if the batteries are installed correctly. The list of possible defects is given below:

1. All symbols are lit at the same time. This happens when the scales are switched on and indicates that the device is in good working order.
2. Err – : the maximum load of the scales has been exceeded.
3. Lo – : the battery is re-energize. Install new batteries.
4. Errh – : the maximum possible for these scales body fat content is exceeded.
5. Errl – : the body fat content is too low.
6. [– : leave the scales and wait a few seconds. Repeat the weighing again.

REFERENCE INFORMATION

The resulting values are based on standard method of measurement. All these indicators are individual and depend on the way of life as well as on variety of indicators. You can see the list of conditions affecting the result obtained in the course of measurements below:

- a) After taking a shower or a bath the indicators can be underestimated.
- b) After eating the indicators can be overestimated.
- c) The indicators can be incorrect during the course of the illness.
- d) You can see values for estimating the amount of body fat:

Age	Fat % Female	Hydration % Female	Fat % Male	Hydration % Male	Status	Key line indication
Below 30	4.0-16.0%	66.0-57.8%	4.0-11.0%	66.0-61.2%	too lean	i
	16.1-20.5%	57.7-54.7%	11.1-15.5%	61.1-58.1%	lean	
	20.6-25.0%	54.6-51.6%	15.6-20.0%	58.0-55.0%	normal	
	25.1-30.5%	51.5-47.8%	20.1-24.5%	54.9-51.9%	high	
	30.6-60.0%	47.7-27.5%	24.6-60.0%	51.8-27.5%	very high	
Above 30	4.0-20.0%	66.0-55.0%	4.0-15.0%	66.0-58.4%	too lean	i
	20.1-25.0%	54.9-51.6%	15.1-19.5%	58.3-55.3%	lean	
	25.1-30.0%	51.5-48.1%	19.6-24.0%	55.2-52.3%	normal	
	30.1-35.0%	48.0-44.7%	24.1-28.5%	52.2-49.2%	high	
	35.1-60.0%	44.6-27.5%	28.6-60.0%	49.1-27.5%	very high	

PLEASE READ FOLLOWING NOTES CAREFULLY FOR PROPER UNDERSTANDING OF THE RESULTS

1. The data obtained during the measurements can be used only as a reference information. Do not use these values for medical purposes. For more detailed information on the recommended diet, etc., contact your doctor.

2. Electronic scales with the function of analysis of the percentage of body fat, hydration and muscle an bone mass can not be used:
- Children under 10, and adults over 80.
 - People with symptoms of fever, edema, osteoporosis.
 - People on dialysis treatment.
 - People taking cardiovascular medicine.
 - Being pregnant.
 - Athletes who undertake more than 10 hours of intensive sport activity.
 - Those with heart pulse rate under 60 beats / minute.

SPECIFICATIONS

Capacity: 150 kg
Division: 100 g
Lifetime Not more than 10000 measurements or
not more than 5 years from the date
of purchase by the buyer.

The head office of the manufacturer : «UNIT Handelsgesellschaft G.m.b.H.»

131, Gersthofer Str., A-1180, Vienna, Austria

Exported by: «UNIT Electronics (H. K.) Limited»

3905 Two Exchange Square, 8 Connaught Place, Central, Hong Kong

Made in People's Republic of China (P.R.C.)

This appliance is to be stored and used under the following conditions: temperature – not less than +5°C; relative humidity - not higher than 80%; aggressive impurities in the air: not allowable. As the Manufacturer is constantly working on the improvement of its products, the design and technical specifications are subject to change without prior notice.

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